

Hospital Mask

NON-VENTED FULL FACE MASK (DISPOSABLE)

Quick-fitting, high-performance

Single-patient, short-term hospital use


- Designed for noninvasive ventilation therapy in the hospital environment.
- Suitable for critical care—effective at high therapy pressures (up to 40 cm H₂O).

Quick and easy to fit and release

- Straight-from-the-box solution; one-piece mask, pre-assembled headgear.
- Quick-release headgear clips—enable quick and easy mask fitting and removal.

Reliable quality and performance


- ResMed's premium dual-wall technology—in a disposable mask.
- High quality—tested to ResMed's rigorous standards.



TECHNOLOGY HIGHLIGHT

Dual-wall cushion technology in a disposable mask

- The thick inner membrane provides support and stability.
- The thin outer membrane increases comfort on contact with the skin.
- Together, the two membranes of the dual-wall cushion offer a gentle, yet stable seal.



Dual-wall cushion technology combines stability (via the thick inner membrane) and comfort (via the thin outer membrane)

PROFESSIONAL BENEFITS

Quick-fitting

- Pre-assembled with headgear
- Designed for minimal adjustment
- Easy-reference fitting template (located on packaging)

Easy-to-use

- One-piece mask frame with no loose parts
- Set-and-forget headgear
- Rotating elbow allows the connecting tube to be positioned in any direction

Efficient inventory management

- Wide fit range with only three sizes

Efficient therapy management

- Optimal seal and minimal leak for effective noninvasive ventilation
- Easy-access port for therapy monitoring and management

Ideal for emergency or critical care

- Quick-release headgear clips save time and allow rapid mask removal if required
- Effective at high therapy pressures (up to 40 cm H₂O)

PATIENT BENEFITS

Increased comfort


- Less pressure on the nasal bridge and reduced chances of skin irritation—the dual-wall cushion is gentle on the face
- Pressure evenly distributed across the head—comfortable headgear design

Convenience

- Quick-release headgear clips allow easy removal of the mask from the headgear without disturbing the patient



Comfortable, 3-strap, set-and-forget headgear: easy and quick to fit; distributes pressure evenly across the head; skin-friendly materials enhance comfort.



Easy to select the right size without opening the bag: the fitting template is integrated into the colour-coded packaging.

PRODUCT CODES

Hospital mask
Disposable non-vented full face mask

ANZ, Asia Pacific

| | |
|----------------------------|--------------|
| Small | 60735 |
| Medium | 60736 |
| Large | 60737 |
| Headgear with clips (10pk) | 60711 |